

AIRVO 2



The Airvo 2 generates high flows of warmed and humidified gases, delivered to the patient through a variety of nasal, tracheostomy and mask interfaces.

Intended Use: The Airvo 2 humidifier is for the treatment of spontaneously breathing patients who would benefit from receiving high flow warmed and humidified respiratory gases. This includes patients who have had upper airways bypassed. The flow may be set from 2Lpm to 60Lpm on the patient interface. The Airvo 2 is for patients in hospitals and long-term care facilities.

UNIT SPECIFICATIONS:

- optional O2 Source, 53 concentrator, OI concentrator or 70lpm flow meeter system.