

# The Vest System



The Hil-Rom Vest System is a unique therapy system designed to help mobilize pulmonary secretions through high frequency chest wall oscillation. The system consists of an inflatable vest or wrap connected by tubes to an air pulse generator. The generator rapidly inflates and deflates the vest, compressing and releasing the chest wall. This technology is called high-frequency chest wall oscillation (HFCWO).

## UNIT SPECIFICATIONS:

- Variable frequency settings: 1-20 HZ
- Variable pressure settings: 1-10 lbs

## OPTIONAL ACCESSORIES:

- X-Small: 19" to 23"
- Small: 23" to 33"
- Medium: 33" to 43"
- Large: 43" to 53"
- X-Large: 53" to 67"
- XX-Large: >67"